gentec-e)

HP125A-15KW-HD-IMP-D0

High power detector for laser power measurement up to 15000 W.



KEY FEATURES

HIGH POWER HANDLING

Handles up to 15 kW of continuous power with our standard models. Custom models available for higher powers (See SUPER HP).

STABLE READING

Less sensitive to variations in water cooling temperature than any other high power water-cooled meter on the market

LARGE APERTURE

Our standard HP models (4KW, 12KW and 15KW) have a very large effective aperture up to 125 mm \emptyset to accomodate even the largest laser beams. Larger apertures with various shapes are available upon request (See SUPER HP).

AVAILABLE WITH YAG AND CO2 CALIBRATION

All HP Models can be calibrated at YAG and $\rm CO_2$ wavelengths with a calibration uncertainty of +/- 5%

DIRECT USB CONNECTION TO A PC

Each head comes with both a DB-15 connector (for use with a Gentec-EO monitor) and a USB2.0 output for direct connection to a PC.

SPECIFICATIONS

MEASUREMENT CAPABILITIES

HE BOREMENT ON MEETIES	
Maximum average power (continuous)	15000 W
Minimum average power ¹	500 W
Noise equivalent power ²	15 W
Spectral range	0.193 - 20 µm
Typical rise time	15 sec
Typical power sensitivity	0.13 mV/W
Power calibration uncertainty	±5 %
Repeatability	±2 %
Back reflections	~ 15 %
Linearity with power	±2 %
Linearity vs beam diameter	±1 %
Linearity vs beam position ³	±1.0 %

1. For lower powers, call your Gentec-EO representative.

2. Nominal value. Actual value depends on electrical noise in the measurement system.

3. For a beam size of 20% of the aperture area, moved across 80% of the aperture area.

WATER REQUIREMENTS

Required cooling flow ¹	(8 - 10) LPM < ±1 LPM/min
Temperature range	15 – 25 °C
Rate of temperature change	< ±3°C/min
Maximum water pressure	413 kPa (60 psi)
1 Contact Gentec-EO for clean dejonized water cooling module option	

1. Contact Gentec-EO for clean deionized water cooling module option

DAMAGE THRESHOLDS

Maximum average power density¹

CONTROLLER AND GUI SPECIFICATIONS	
Data display	Real time, scope, needle, averaging, histogram and statistics
Analog output ¹	0-2 Volts
Serial commands via	USB
External power supply ²	Through USB or Gentec-EO displays & PC interfaces
Display type	None
 12 V maximum output signal available upon request. 2. A USB power adaptor will be necessary if the HP is used with a DB-15 extension cable. 	
PHYSICAL CHARACTERISTICS	
Aperture width	125 mm
Aperture height	125 mm
Absorber	HD
Dimensions	153H x 153W x 70D mm
Weight	5 kg
ORDERING INFORMATIONS	
HP125A-15KW-HD-IMP-D0	2022631
HP125A-15KW-HD-IMP-BLU-D0	TBD

INTERESTED IN THIS PRODUCT?

GET A QUOTE

Find your local sales representative at gentec-eo.com/contact-us